


ARL4A Antibody (R34114)

Catalog No.	Formulation	Size
R34114-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Pig, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	10124
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:1000
Limitations	This ARL4A antibody is available for research use only.



Western blot testing of human testis lysate with ARL4A antibody at 2ug/ml.
Predicted/observed molecular weight: ~23kDa.

Description

Additional name(s) for this target protein: ADP-ribosylation factor-like 4

Application Notes

Optimal dilution of the ARL4A antibody should be determined by the researcher.

Immunogen

Amino acids GNGLSDQTSILSN were used as the immunogen for this ARL4A antibody.

Storage

Aliquot and store the ARL4A antibody at -20oC.