


ARL6IP5 Antibody (R34905)


Catalog No.	Formulation	Size
R34905-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)


Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	10550
Localization	Cytoplasmic
Applications	Western Blot : 0.1-0.3ug/ml IHC (FFPE) : 3.75ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This ARL6IP5 antibody is available for research use only.



Western blot testing of human brain lysate with ARL6IP5 antibody at 0.1ug/ml. Predicted molecular weight: ~22 kDa.



IHC testing of FFPE human small intestine tissue with ARL6IP5 antibody at 3.75ug/ml. Required HIER: steamed antigen retrieval with pH6 citrate buffer; AP-staining.



IHC testing of FFPE human tonsil tissue with ARL6IP5 antibody at 3.75ug/ml. Required HIER: steamed antigen retrieval with pH6 citrate buffer; AP-staining.

Description

Additional name(s) for this target protein: ADP-ribosylation-like factor 6 interacting protein 5, GTRAP3-18, JWA

Application Notes

Optimal dilution of the ARL6IP5 antibody should be determined by the researcher.

Immunogen

Amino acids NRLTDYISKVKE were used as the immunogen for this ARL6IP5 antibody.

Storage

Aliquot and store the ARL6IP5 antibody at -20oC.