


ARMET Antibody (R34428)

Catalog No.	Formulation	Size
R34428-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	7873
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This ARMET antibody is available for research use only.



Western blot testing of human pancreas lysate with ARMET antibody at 0.3ug/ml.
Predicted molecular weight: ~21kDa.

Description

Additional name(s) for this target protein: Arginine-rich, mutated in early stage tumors

Application Notes

Optimal dilution of the ARMET antibody should be determined by the researcher.

Immunogen

Amino acids KFCREARGKENR were used as the immunogen for this ARMET antibody.

Storage

Aliquot and store the ARMET antibody at -20oC.