


Aryl Hydrocarbon Receptor Antibody (R34460)

| Catalog No. | Formulation | Size |
|--------------|---|--------|
| R34460-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

| | |
|--------------------|---|
| Availability | 1-3 business days |
| Species Reactivity | Human |
| Format | Antigen affinity purified |
| Clonality | Polyclonal (goat origin) |
| Isotype | Goat Ig |
| Purity | Antigen affinity |
| Gene ID | 196 |
| Applications | Western Blot : 0.2-0.5ug/ml ELISA (peptide) LOD : 1:128000 |
| Limitations | This Aryl Hydrocarbon Receptor antibody is available for research use only. |



Western blot testing of HeLa lysate with Aryl Hydrocarbon Receptor antibody at 0.2ug/ml. Predicted molecular weight ~ 95kDa.

Description

Additional name(s) for this target protein: AH-receptor; AHR

Application Notes

Optimal dilution of the Aryl Hydrocarbon Receptor antibody should be determined by the researcher.

Immunogen

Amino acids KHGLNPQSAITPQT were used as the immunogen for this Aryl Hydrocarbon Receptor antibody.

Storage

Aliquot and store the Aryl Hydrocarbon Receptor antibody at -20oC.