


ATG4A Antibody (R36152)

Catalog No.	Formulation	Size
R36152-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	115201
Applications	Immunohistochemistry (FFPE) : 2-4ug/ml ELISA (peptide) LOD : 1:1000
Limitations	This ATG4A antibody is available for research use only.



IHC testing of FFPE human liver with ATG4A antibody at 2ug/ml. HIER: steamed with pH6 citrate buffer, HRP-staining.

Description

Additional name(s) for this target protein: Autophagy related 4A

Application Notes

Optimal dilution of the ATG4A antibody should be determined by the researcher.

Immunogen

Amino acids TEENGTVNDQTFHC were used as the immunogen for this ATG4A antibody.

Storage

Aliquot and store the ATG4A antibody at -20oC.