


ATG4C Antibody / APG4C (R33726)

Catalog No.	Formulation	Size
R33726-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Cow
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	84938
Applications	Western Blot : 2-6ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This ATG4C antibody is available for research use only.



Western blot testing of human skeletal muscle lysate with ATG4C antibody at 2ug/ml.
Predicted molecular weight: ~52kDa.

Description

Additional name(s) for this target protein: ATG4 autophagy related 4 homolog C

Application Notes

Optimal dilution of the ATG4C antibody should be determined by the researcher.

Immunogen

Amino acids EDEKKQLKRFSTEE were used as the immunogen for this ATG4C antibody.

Storage

Aliquot and store the ATG4C antibody at -20oC.