


ATP5C1 Antibody (R36220)

Catalog No.	Formulation	Size
R36220-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Rat
Predicted Reactivity	Cow, Dog, Mouse
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	509
Applications	Western Blot : 0.01-0.03ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This ATP5C1 antibody is available for research use only.



Western blot testing of rat heart lysate with ATP5C1 antibody at 0.01ug/ml.
Predicted/observed molecular weight: ~33kDa.

Description

Additional name(s) for this target protein: ATP synthase, H⁺ transporting, mitochondrial F1 complex, gamma polypeptide 1; ATPG

Application Notes

Optimal dilution of the ATP5C1 antibody should be determined by the researcher.

Immunogen

Amino acids TLTFNRTRQAVITKE were used as the immunogen for this ATP5C1 antibody.

Storage

Aliquot and store the ATP5C1 antibody at -20oC.