


ATRIP Antibody / ATR-interacting protein [clone 26A07] (RQ8899)

Catalog No.	Formulation	Size
RQ8899	Antibody in PBS with 0.02% sodium azide, 50% glycerol and 0.4-0.5mg/ml BSA	100 ul


Recombinant **RABBIT MONOCLONAL**

[Bulk quote request](#)

Availability	1-3 days
Species Reactivity	Human
Format	Purified
Clonality	Recombinant Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	26A07
Purity	Affinity chromatography
UniProt	Q8WXE1
Localization	Nuclear
Applications	Western Blot : 1:500 Immunohistochemistry (FFPE) : 1:50
Limitations	This ATRIP antibody is available for research use only.



IHC staining of FFPE human lung cancer tissue with ATRIP antibody, HRP-secondary and DAB substrate. HIER: boil tissue sections in pH8 EDTA for 20 min and allow to cool before testing.



Western blot testing of human 1) PC-3, 2) HEL and 3) U-251 cell lysate with ATRIP antibody. Predicted molecular weight ~86 kDa.

