


AXIN1 Antibody (R34056)

Catalog No.	Formulation	Size
R34056-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	8312
Applications	Immunohistochemistry (FFPE) : 3-5ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This AXIN1 antibody is available for research use only.



IHC testing of FFPE human adrenal cortex with AXIN1 antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

Description

Additional name(s) for this target protein: Axis inhibition protein 1

Application Notes

Optimal dilution of the AXIN1 antibody should be determined by the researcher.

Immunogen

Amino acids FEEKIIGKVEKVD were used as the immunogen for this AXIN1 antibody.

Storage

Aliquot and store the AXIN1 antibody at -20oC.