


## B7-H4 Antibody (R35936)

Catalog No.	Formulation	Size
R35936-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	79679
Applications	Western Blot : 0.5-1.5ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This B7-H4 antibody is available for research use only.



Western blot testing of human pancreas lysate with B7-H4 antibody at 0.5ug/ml.  
Predicted molecular weight ~31kDa.

## Description

Additional name(s) for this target protein: B7 family member, H4; VTCN1; V-set domain containing T cell activation inhibitor 1

## Application Notes

Optimal dilution of the B7-H4 antibody should be determined by the researcher.

## Immunogen

Amino acids SKGKGNNANLEYK were used as the immunogen for this B7-H4 antibody.

## **Storage**

Aliquot and store the B7-H4 antibody at -20oC.