


BAG3 Antibody [clone ACBO-2] (RQ4815)

Catalog No.	Formulation	Size
RQ4815	Antibody in PBS with 0.02% sodium azide, 50% glycerol and 0.4-0.5mg/ml BSA	100 ul

[Bulk quote request](#)

Availability	1-2 weeks
Species Reactivity	Human
Format	Purified
Clonality	Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	ACBO-2
Purity	Affinity purified
UniProt	O95817
Applications	Western Blot : 1:500-1:2000
Limitations	This BAG3 antibody is available for research use only.



Western blot testing of human K562 cell lysate with BAG3 antibody. Predicted molecular weight ~62 kDa but can be observed at 75~80 kDa.

Description

BAG family molecular chaperone regulator 3 inhibits the chaperone activity of HSP70/HSC70 by promoting substrate release. Has anti-apoptotic activity.

Application Notes

Optimal dilution of the BAG3 antibody should be determined by the researcher.

Immunogen

A synthetic peptide specific to human BAG3 was used as the immunogen for the BAG3 antibody.

Storage

Store the BAG3 antibody at -20oC.