


BAG5 Antibody (R33549)


Catalog No.	Formulation	Size
R33549-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat, Dog, Cow
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	9529
Localization	Cytoplasmic, nuclear
Applications	Western Blot : 1-3ug/ml IHC (FFPE) : 3.75ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This BAG5 antibody is available for research use only.



Western blot testing of nuclear HeLa lysate with BAG5 antibody at 2ug/ml. Predicted molecular weight: ~56/51 kDa (isoforms a/b).



IHC testing of FFPE human testis tissue with BAG5 antibody at 3.75ug/ml. Required HIER: steamed antigen retrieval with pH6 citrate buffer; AP-staining.

Description

Additional name(s) for this target protein: BCL2-associated athanogene 5

Application Notes

Optimal dilution of the BAG5 antibody should be determined by the researcher.

Immunogen

Amino acids DMGNQHPSISRLQ were used as the immunogen for this BAG5 antibody.

Storage

Aliquot and store the BAG5 antibody at -20oC.