


BAG5 Antibody (R33579)

Catalog No.	Formulation	Size
R33579-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	9529
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This BAG5 antibody is available for research use only.



Western blot testing of Jurkat lysate with BAG5 antibody at 1ug/ml. Predicted molecular weight: ~56/51kDa (isoforms a/b).

Description

Additional name(s) for this target protein: BCL2-associated athanogene 5

Application Notes

Optimal dilution of the BAG5 antibody should be determined by the researcher.

1. This BAG5 antibody is specific for isoforms a and b of the protein.

Immunogen

Amino acids DPQGEEKCKAARKQ were used as the immunogen for this BAG5 antibody.

Storage

Aliquot and store the BAG5 antibody at -20oC.