


Bcl2a1a Antibody (R35324)

Catalog No.	Formulation	Size
R35324-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
UniProt	Q07440
Gene ID	12044
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This Bcl2a1a antibody is available for research use only.



Western blot testing of rat spleen lysate with Bcl2a1a antibody at 0.1ug/ml. Predicted molecular weight: ~20kDa but may be observed at up to ~37kDa.

Description

Additional name(s) for this target protein: Bcl-2-related protein A1

Application Notes

Optimal dilution of the Bcl2a1a antibody should be determined by the researcher.

1. This antibody may cross react with the highly homologous Bcl2a1b and Bcl2a1d.

Immunogen

Amino acids QKEVEKNLKPYLDD were used as the immunogen for this Bcl2a1a antibody.(1)

Storage

Aliquot and store the Bcl2a1a antibody at -20oC.