


## Beclin Antibody (F46280)


| Catalog No.   | Formulation                                | Size    |
|---------------|--------------------------------------------|---------|
| F46280-0.4ML  | In 1X PBS, pH 7.4, with 0.09% sodium azide | 0.4 ml  |
| F46280-0.08ML | In 1X PBS, pH 7.4, with 0.09% sodium azide | 0.08 ml |

[Bulk quote request](#)

|                             |                                                                                    |
|-----------------------------|------------------------------------------------------------------------------------|
| <b>Availability</b>         | 1-3 business days                                                                  |
| <b>Species Reactivity</b>   | Human, Mouse                                                                       |
| <b>Predicted Reactivity</b> | Bovine, Pig, Rat                                                                   |
| <b>Format</b>               | Purified                                                                           |
| <b>Clonality</b>            | Polyclonal (rabbit origin)                                                         |
| <b>Isotype</b>              | Rabbit Ig                                                                          |
| <b>Purity</b>               | Purified                                                                           |
| <b>UniProt</b>              | Q14457                                                                             |
| <b>Applications</b>         | Western Blot : 1:1000<br>IHC (Paraffin) : 1:50-1:100<br>Flow Cytometry : 1:10-1:50 |
| <b>Limitations</b>          | This Beclin antibody is available for research use only.                           |




Western blot testing of Beclin antibody and (1) A2058, (2) HeLa, (3) mouse brain, (4) Y79, and (5) HL-60 lysate



Western blot analysis of Beclin in 293 lysate transiently transfected with the BECN1 gene (2ug/lane).



IHC analysis of FFPE human breast carcinoma tissue stained with the Beclin antibody



Flow cytometric analysis of Y79 cells using Beclin antibody (bottom histogram) compared to a negative control (top histogram). FITC-conjugated goat-anti-rabbit secondary Ab was used for the analysis.

## Description

Macroautophagy is the major inducible pathway for the general turnover of cytoplasmic constituents in eukaryotic cells, it is also responsible for the degradation of active cytoplasmic enzymes and organelles during nutrient starvation. Macroautophagy involves the formation of double-membrane bound autophagosomes which enclose the cytoplasmic constituent targeted for degradation in a membrane bound structure, which then fuse with the lysosome (or vacuole) releasing a single-membrane bound autophagic bodies which are then degraded within the lysosome (or vacuole). Beclin 1 plays a role in two fundamentally important cell biological pathways: autophagy and apoptosis. Beclin 1 is thought to function as a VPS and autophagy protein as part of a complex with Class III PI3 kinase, Vps34.

## Application Notes

Titration of the Beclin antibody may be required due to differences in protocols and secondary/substrate sensitivity.

## Immunogen

A portion of amino acids 210-239 from the human protein was used as the immunogen for this Beclin antibody.

## Storage

Aliquot the Beclin antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.