


Beta-Arrestin 2 Antibody (R34167)


Catalog No.	Formulation	Size
R34167-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
UniProt	P32121
Gene ID	409
Applications	Western Blot : 0.3-1ug/ml Immunohistochemistry (FFPE) : 4-6ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This beta-Arrestin 2 antibody is available for research use only.



Western blot testing of human brain lysate with beta-Arrestin 2 antibody at 0.3ug/ml. Predicted molecular weight: 44~47 kDa (isoforms 1-5).



IHC staining of FFPE human cerebellum with beta-Arrestin 2 antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, HRP-staining.

Description

Additional name(s) for this target protein: ARB2, ARR2; ARRB2

Application Notes

Optimal dilution of the beta-Arrestin 2 antibody should be determined by the researcher.

Immunogen

Amino acids HDHIPLPRPQS were used as the immunogen for this beta-Arrestin 2 antibody. This sequence is common to isoforms 1-5 of P32121 (UniProt).

Storage

Aliquot and store the beta-Arrestin 2 antibody at -20oC.