


## Beta-catenin Antibody (R35308)

Catalog No.	Formulation	Size
R35308-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Pig, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	1499
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This Beta-catenin antibody is available for research use only.



Western blot testing of HeLa lysate with Beta-catenin antibody at 0.3ug/ml. Predicted molecular weight ~85 kDa, but routinely observed at 90-95 kDa.

## Description

Additional name(s) for this target protein: Catenin beta-1; CTNNB1

## Application Notes

Optimal dilution of the Beta-catenin antibody should be determined by the researcher.

## Immunogen

Amino acids QIPSTQFDAAHPTN were used as the immunogen for this Beta-catenin antibody.

## Storage

Aliquot and store the Beta-catenin antibody at -20oC.