


## BHLHE22 Antibody (R36464)

| Catalog No.  | Formulation                                                                   | Size   |
|--------------|-------------------------------------------------------------------------------|--------|
| R36464-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

|                    |                                                           |
|--------------------|-----------------------------------------------------------|
| Availability       | 1-3 business days                                         |
| Species Reactivity | Human, Pig                                                |
| Format             | Antigen affinity purified                                 |
| Clonality          | Polyclonal (goat origin)                                  |
| Isotype            | Goat Ig                                                   |
| Purity             | Antigen affinity                                          |
| UniProt            | Q8NFJ8                                                    |
| Gene ID            | 27319                                                     |
| Applications       | Western Blot : 1-3ug/ml<br>ELISA (peptide) LOD : 1:32000  |
| Limitations        | This BHLHE22 antibody is available for research use only. |



Western blot testing of human cerebellum lysate with BHLHE22 antibody at 1ug/ml.  
Predicted molecular weight ~37 kDa.

## Description

Additional name(s) for this target protein: Basic helix-loop-helix family, member e22; bHLHb5; Trinucleotide repeat-containing gene 20 protein; TNRC20

## Application Notes

Optimal dilution of the BHLHE22 antibody should be determined by the researcher.

## Immunogen

Amino acids KYGESASRGSSVAES were used as the immunogen for this BHLHE22 antibody.

## **Storage**

Aliquot and store the BHLHE22 antibody at -20oC.