


BLNK Antibody (R34385)


Catalog No.	Formulation	Size
R34385-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Pig, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	29760
Applications	Western Blot : 4-8ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This BLNK antibody is available for research use only.



Western blot testing of Daudi lysate with BLNK antibody at 4ug/ml. Predicted molecular weight: ~51kDa but can be observed at 65-70kDa.



Western blot of HEK293 lysate overexpressing BLNK probed with BLNK antibody (mock transfection in lane 1). Predicted molecular weight: ~51kDa but can be observed at 65-70kDa.

Description

Additional name(s) for this target protein: B-cell linker, BASH, SLP-65, Src homology 2 domain-containing leukocyte protein of 65 kDa

Application Notes

Optimal dilution of the BLNK antibody should be determined by the researcher.

Immunogen

Amino acids KDSTRLKYAVKVS were used as the immunogen for this BLNK antibody.

Storage

Aliquot and store the BLNK antibody at -20oC.