


BMP15 Antibody [clone BMP15/4263] (V8690)

Catalog No.	Formulation	Size
V8690-100UG	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide	100 ug
V8690-20UG	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide	20 ug
V8690SAF-100UG	1 mg/ml in 1X PBS; BSA free, sodium azide free	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Clonality	Monoclonal (mouse origin)
Isotype	Mouse IgG1
Clone Name	BMP15/4263
Purity	Protein G affinity chromatography
UniProt	O95972
Localization	Cytoplasmic, secreted
Applications	ELISA : order Ab without BSA for coating Western Blot : 1-2ug/ml
Limitations	This BMP15 antibody is available for research use only.



SDS-PAGE analysis of purified, BSA-free BMP15 antibody as confirmation of integrity and purity.

Description

BMP15, also known as GDF9B, plays a crucial role in the regulation of fertility.

Application Notes

Optimal dilution of the BMP15 antibody should be determined by the researcher.

Immunogen

Amino acids SAEVTASSSKHSGPENNC from the C-terminus of the human protein were used as the immunogen for the BMP15 antibody.

Storage

Store the BMP15 antibody at 2-8oC (with azide) or aliquot and store at -20oC or colder (without azide).