


## BMP7 Antibody (R36451)

| Catalog No.  | Formulation                                                                   | Size   |
|--------------|-------------------------------------------------------------------------------|--------|
| R36451-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

|                    |                                                         |
|--------------------|---------------------------------------------------------|
| Availability       | 1-3 business days                                       |
| Species Reactivity | Human, Mouse, Rat, Dog, Pig, Cow                        |
| Format             | Antigen affinity purified                               |
| Clonality          | Polyclonal (goat origin)                                |
| Isotype            | Goat Ig                                                 |
| Purity             | Antigen affinity                                        |
| UniProt            | P18075                                                  |
| Gene ID            | 655                                                     |
| Applications       | Western Blot : 1-3ug/ml<br>ELISA (peptide) LOD : 1:2000 |
| Limitations        | This BMP7 antibody is available for research use only.  |



Western blot testing of pig kidney lysate with BMP7 antibody at 1ug/ml. Predicted molecular weight: ~49 kDa.

## Description

Additional name(s) for this target protein: Bone morphogenetic protein 7

## Application Notes

Optimal dilution of the BMP7 antibody should be determined by the researcher.

## Immunogen

Amino acids 120-131 (SLQDSHFLTDAD) were used as the immunogen for this BMP7 antibody.

## Storage

Aliquot and store the BMP7 antibody at -20oC.