


BMP8A Antibody (F47545)

Catalog No.	Formulation	Size
F47545-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F47545-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	Q7Z5Y6
Applications	Western Blot : 1:1000
Limitations	This BMP8A antibody is available for research use only.



Description

Induces cartilage and bone formation. May be the osteoinductive factor responsible for the phenomenon of epithelial osteogenesis. Plays a role in calcium regulation and bone homeostasis (By similarity). [UniProt]

Application Notes

Titration of the BMP8A antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 271-298 from the human protein was used as the immunogen for this BMP8A antibody.

Storage

Aliquot the BMP8A antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.