


BNIP1 Antibody (F47254)

Catalog No.	Formulation	Size
F47254-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F47254-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	Q12981
Applications	Western Blot : 1:1000
Limitations	This BNIP1 antibody is available for research use only.



Description

This gene is a member of the BCL2/adenovirus E1B 19 kd-interacting protein (BNIP) family. It interacts with the E1B 19 kDa protein which is responsible for the protection of virally-induced cell death, as well as E1B 19 kDa-like sequences of BCL2, also an apoptotic protector. Alternative splicing of this gene results in four protein products with identical N- and C-termini.

Application Notes

Titration of the BNIP1 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 26-55 from the human protein was used as the immunogen for this BNIP1 antibody.

Storage

Aliquot the BNIP1 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.