


BNIP1 Antibody (R33537)

Catalog No.	Formulation	Size
R33537-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	662
Applications	Western Blot : 0.5-1.5ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This BNIP1 antibody is available for research use only.



Western blot testing of human kidney lysate with BNIP1 antibody at 0.5ug/ml. Predicted molecular weight: ~26kDa. Both bands are both blocked by the immunizing peptide.

Description

Additional name(s) for this target protein: BCL2/adenovirus E1B 19kDa interacting protein 1

Application Notes

Optimal dilution of the BNIP1 antibody should be determined by the researcher.

Immunogen

Amino acids DLLRQRKTTKESLAQ were used as the immunogen for this BNIP1 antibody.

Storage

Aliquot and store the BNIP1 antibody at -20oC.