


BNIP3L Antibody / Nix (F53746)

Catalog No.	Formulation	Size
F53746-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F53746-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Format	Purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Purified
UniProt	O60238
Applications	IHC (Paraffin) : 1:50-1:100
Limitations	This BNIP3L antibody is available for research use only.



IHC analysis of FFPE human hepatocarcinoma tissue stained with the BNIP3L antibody

Description

BNIP3L is a member of the BCL2/adenovirus E1B 19 kd-interacting protein (BNIP) family. It interacts with the E1B 19 kDa protein which is responsible for the protection of virally-induced cell death, as well as E1B 19 kDa-like sequences of BCL2, also an apoptotic protector. The protein encoded by this gene is a functional homolog of BNIP3, a proapoptotic protein. This protein may function simultaneously with BNIP3 and may play a role in tumor suppression.

Application Notes

Titration of the BNIP3L antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 122-152 from the human protein (BH3 domain) was used as the immunogen for this BNIP3L antibody.

Storage

Aliquot the BNIP3L antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.