


Bradykinin Antibody (R35131)

Catalog No.	Formulation	Size
R35131-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3827
Applications	Western Blot : 0.03-0.1ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This Bradykinin antibody is available for research use only.



Western blot testing of albumin-depleted plasma lysate with Bradykinin antibody at 0.03ug/ml. Predicted/observed molecular weight ~72kDa.

Description

Additional name(s) for this target protein: BDK, KNG1, kininogen 1; KNG1

Application Notes

Optimal dilution of the Bradykinin antibody should be determined by the researcher.

Immunogen

Amino acids PRDIPTNSPELE were used as the immunogen for this Bradykinin antibody.

Storage

Aliquot and store the Bradykinin antibody at -20oC.