


c-Met Antibody [clone 8C10-F6-C11] (F54057)

Catalog No.	Formulation	Size
F54057-0.1ML	In PBS with 50% glycerol and 0.03% ProClin 300	0.1 ml

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Clonality	Monoclonal (mouse origin)
Isotype	Mouse IgG2a
Clone Name	8C10-F6-C11
Purity	Protein G affinity
UniProt	P08581
Applications	Western Blot : 1:1000 (transfected lysate)
Limitations	This c-Met antibody is available for research use only.



Western blot testing of c-Met A) transfected (partial protein) and (B) untransfected CHO-K1 cell lysate with c-Met antibody at 1:1000. Predicted molecular weight of the transfected partial: ~38 kDa.

Description

Receptor for hepatocyte growth factor and scatter factor. Has a tyrosine-protein kinase activity. Functions in cell proliferation, scattering, morphogenesis and survival.

Application Notes

The stated application concentrations are suggested starting points. Titration of the c-Met antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A human recombinant protein was used as the immunogen for this c-Met antibody.

Storage

Store the c-Met antibody at -20oC.