


C18orf32 Antibody (F53867)

Catalog No.	Formulation	Size
F53867-0.05ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.05 ml

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	Q8TCD1
Applications	Western Blot : 1:1000-2000
Limitations	This C18orf32 antibody is available for research use only.



Western blot testing of human 1) human brain and 2) LNCaP cell lysate with C18orf32 antibody at 1:2000. Predicted molecular weight: 9 kDa.

Description

May activate the NF-kappa-B signaling pathway.

Application Notes

Titration of the C18orf32 antibody may be required due to differences in protocols and substrate sensitivity.

Immunogen

A portion of amino acids 26-58 from the human protein was used as the immunogen for the C18orf32 antibody.

Storage

Aliquot the C18orf32 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.

