


C1D Antibody (R35942)

Catalog No.	Formulation	Size
R35942-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	10438
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This C1D antibody is available for research use only.



Western blot testing of human skeletal muscle lysate with C1D antibody at 0.1ug/ml.
Predicted molecular weight: ~16kDa.

Description

Additional name(s) for this target protein: Nuclear nucleic acid-binding protein C1D; Small unique nuclear receptor corepressor (SUNCOR)

Application Notes

Optimal dilution of the C1D antibody should be determined by the researcher.

Immunogen

Amino acids SKNASKVANKGKSKS were used as the immunogen for this C1D antibody.

Storage

Aliquot and store the C1D antibody at -20oC.