

C5 Antibody (F51211)


Catalog No.	Formulation	Size
F51211-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F51211-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)


Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	P01031
Applications	Western Blot : 1:1000 Immunohistochemistry (Paraffin) : 1:50-1:100 Flow Cytometry : 1:10-1:50
Limitations	This C5 antibody is available for research use only.



Western blot analysis of C5 antibody and HL-60 lysate



IHC analysis of FFPE human breast carcinoma with C5 antibody



Flow cytometric analysis of MDA-231 cells using C5 antibody (green) compared to a