


C8B Antibody / Complement C8B (F47544)

Catalog No.	Formulation	Size
F47544-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F47544-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	P07358
Applications	Western Blot : 1:1000
Limitations	This C8B antibody is available for research use only.



Description

C8 beta is one of the three subunits that comprise the component 8 (C8) of the complement system. C8 participates in the formation of Membrane Attack Complex that results in the lysis of cells. Patients with C8B deficiency are prone to bacteria infection.

Application Notes

Titration of the C8B antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 393-422 from the human protein was used as the immunogen for this C8B antibody.

Storage

Aliquot the C8B antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.