


Cadherin 23 Antibody (R36389)

Catalog No.	Formulation	Size
R36389-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	64072
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:1000
Limitations	This Cadherin 23 antibody is available for research use only.



Western blot testing of human amygdala lysate with Cadherin 23 antibody at 1ug/ml.
Predicted/observed molecular weight: ~60kDa.

Description

Additional name(s) for this target protein: CDH23

Application Notes

Optimal dilution of the Cadherin 23 antibody should be determined by the researcher.

Immunogen

Amino acids YNISLYENVTVGTS were used as the immunogen for this Cadherin 23 antibody.

Storage

Aliquot and store the Cadherin 23 antibody at -20oC.