


Calpain 1 Antibody (R34767)

Catalog No.	Formulation	Size
R34767-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Pig
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	823
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This Calpain 1 antibody is available for research use only.



Western blot testing of human kidney lysate with Calpain 1 antibody at 0.3ug/ml.
Predicted molecular weight: ~82kDa.

Description

Additional name(s) for this target protein: CAPN1

Application Notes

Optimal dilution of the Calpain 1 antibody should be determined by the researcher.

Immunogen

Amino acids NASRRARSEQFINLR were used as the immunogen for this Calpain 1 antibody.

Storage

Aliquot and store the Calpain 1 antibody at -20oC.