


Calsenilin Antibody (R35965)

Catalog No.	Formulation	Size
R35965-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	30818
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:1000
Limitations	This Calsenilin antibody is available for research use only.



Western blot of HEK293 lysate overexpressing human KCNIP3-MYC probed with CSEN antibody [1ug/ml] in Lane A and anti-MYC [1/1000] in lane C. Mock-transfected HEK293 probed with CSEN antibody [1ug/ml] in Lane B. Predicted molecular weight: ~29kDa.

Description

Additional name(s) for this target protein: CSEN, Kv channel interacting protein 3, KCNIP3, DREAM

Application Notes

Optimal dilution of the Calsenilin antibody should be determined by the researcher.

1. This Calsenilin antibody is specific for isoform 1.

Immunogen

Amino acids SLLGDLGHTPLS were used as the immunogen for this calsenilin antibody.

Storage

Aliquot and store the Calsenilin antibody at -20oC.