


CAPG Antibody (R34825)


Catalog No.	Formulation	Size
R34825-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	822
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This CAPG antibody is available for research use only.



Western blot testing of U937 [A] and human spleen [B] lysates with CAPG antibody at 0.3ug/ml. Predicted molecular weight: ~38kDa.



Western blot testing of peripheral blood lymphocyte lysate with CAPG antibody at 0.3ug/ml (Lane A: cat # R34825 and Lane B: [cat # 34356](#)). Predicted molecular weight: ~38kDa.

Description

Additional name(s) for this target protein: Capping protein (actin filament), gelsolin-like, Macrophage capping protein, MCP; CAPG

Application Notes

Optimal dilution of the CAPG antibody should be determined by the researcher.

Immunogen

Amino acids NIRATERALNWDS were used as the immunogen for this CAPG antibody.

Storage

Aliquot and store the CAPG antibody at -20oC.