


CASC3 Antibody (F40083)


Catalog No.	Formulation	Size
F40083-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F40083-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat, Bovine, Zebrafish, Xenopus
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	O15234
Applications	IHC (Paraffin) : 1:10-1:50 Western Blot : 1:1000
Limitations	This CASC3 antibody is available for research use only.



IHC analysis of FFPE human breast carcinoma tissue stained with CASC3 antibody



Western blot analysis of CASC3 antibody and NCI-H460 lysate

Description

Component of the dendritic ribonucleoprotein particles (RNPs) in hippocampal neurons. May play a role in mRNA transport.

Application Notes

Titration of the CASC3 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 159-188 from the human protein was used as the immunogen for this CASC3 antibody.

Storage

Aliquot the CASC3 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.