


## Caspase 3 Antibody (R35575)

Catalog No.	Formulation	Size
R35575-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	836
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This Caspase 3 antibody is available for research use only.



Western blot testing of MOLT4 lysate with Caspase 3 antibody at 1ug/ml. The expected ~32kDa band and the additional ~65kDa band are both blocked by the immunizing peptide.

## Description

Additional name(s) for this target protein: CASP3, SCA-1, YAMA, Apopain

## Application Notes

Optimal dilution of the Caspase 3 antibody should be determined by the researcher.

## Immunogen

Amino acids RDVSKEDHSKRS were used as the immunogen for this Caspase 3 antibody.

## Storage

Aliquot and store the Caspase 3 antibody at -20oC.