


CB1 Antibody (R34829)

Catalog No.	Formulation	Size
R34829-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Predicted Reactivity	Cow, Dog, Pig
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	1268
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This CB1 antibody is available for research use only.



Western blot testing of rat brain lysate with CB1 antibody at 0.3ug/ml. Predicted molecular weight ~53kDa, can also be observed at ~70kDa.

Description

Additional name(s) for this target protein: Cannabinoid receptor 1 (brain); CNR1

Application Notes

Optimal dilution of the CB1 antibody should be determined by the researcher.

Immunogen

Amino acids NITEFYNKSLSSFKE were used as the immunogen for this CB1 antibody.

Storage

Aliquot and store the CB1 antibody at -20oC.