


## CCL8 Antibody / MCP2 [clone CCL8/3683] (V9560)


Catalog No.	Formulation	Size
V9560-100UG	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide	100 ug
V9560-20UG	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide	20 ug
V9560SAF-100UG	1 mg/ml in 1X PBS; BSA free, sodium azide free	100 ug

[Bulk quote request](#)

<b>Availability</b>	1-3 business days
<b>Species Reactivity</b>	Human
<b>Format</b>	Purified
<b>Clonality</b>	Monoclonal (mouse origin)
<b>Isotype</b>	Mouse IgG2c, kappa
<b>Clone Name</b>	CCL8/3683
<b>Purity</b>	Protein A/G affinity
<b>UniProt</b>	P80075
<b>Localization</b>	Secreted
<b>Applications</b>	Immunohistochemistry (FFPE) : 1-2ug/ml
<b>Limitations</b>	This CCL8 antibody is available for research use only.



IHC staining of FFPE human adrenal gland with CCL8 antibody (clone CCL8/3683).  
 HIER: boil tissue sections in pH 9 10mM Tris with 1mM EDTA for 20 min and allow to cool before testing.



SDS-PAGE analysis of purified, BSA-free CCL8 antibody (clone CCL8/3683) as confirmation of integrity and purity.

## Description

This anti-microbial gene is one of several chemokine genes clustered on the q-arm of chromosome 17. Chemokines form a superfamily of secreted proteins involved in immunoregulatory and inflammatory processes. The superfamily is divided into four subfamilies based on the arrangement of N-terminal cysteine residues of the mature peptide. This chemokine is a member of the CC subfamily which is characterized by two adjacent cysteine residues. This cytokine displays chemotactic activity for monocytes, lymphocytes, basophils and eosinophils. By recruiting leukocytes to sites of inflammation this cytokine may contribute to tumor-associated leukocyte infiltration and to the antiviral state against HIV infection.

## Application Notes

Optimal dilution of the CCL8 antibody should be determined by the researcher.

## Immunogen

A portion of amino acids 1-99 from the human protein was used as the immunogen for the CCL8 antibody.

## Storage

Aliquot the CCL8 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.