


CCR10 Antibody (RQ4906)


Catalog No.	Formulation	Size
RQ4906	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

Bulk quote request


Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit IgG
Purity	Antigen affinity purified
Buffer	Lyophilized from 1X PBS with 2% Trehalose and 0.025% sodium azide
UniProt	P46092
Applications	Western Blot : 0.5-1ug/ml Immunohistochemistry (FFPE) : 1-2ug/ml
Limitations	This CCR10 antibody is available for research use only.




IHC staining of FFPE human esophageal squamous cancer with CCR10 antibody at 1ug/ml. HIER: boil tissue sections in pH6, 10mM citrate buffer, for 10-20 min and allow to cool before testing.



IHC staining of FFPE human prostate cancer with CCR10 antibody at 1ug/ml. HIER: boil tissue sections in pH6, 10mM citrate buffer, for 10-20 min and allow to cool before testing.



IHC staining of FFPE human tonsil with CCR10 antibody at 1ug/ml. HIER: boil tissue sections in pH6, 10mM citrate buffer, for 10-20 min and allow to cool before testing.



Western blot testing of human 1) HeLa, 2) K562, 3) PC-3, 4) A549 and 5) T-47D lysate with CCR10 antibody at 0.5ug/ml. Predicted molecular weight ~38 kDa.

Description

C-C chemokine receptor type 10 is a protein that in humans is encoded by the CCR10 gene. It is mapped to 17q21.2. Chemokines are a group of small (approximately 8 to 14 kD), mostly basic, structurally related molecules that regulate cell trafficking of various types of leukocytes through interactions with a subset of 7-transmembrane, G protein-coupled receptors. Chemokines also play fundamental roles in the development, homeostasis, and function of the immune system, and they have effects on cells of the central nervous system as well as on endothelial cells involved in angiogenesis or angiostasis. Chemokines are divided into 2 major subfamilies, CXC and CC, based on the arrangement of the first 2 of the 4 conserved cysteine residues; the 2 cysteines are separated by a single amino acid in CXC chemokines and are adjacent in CC chemokines. CCR10 is the receptor for CCL27 (SCYA27; MIM 604833); CCR10-CCL27 interactions are involved in T cell-mediated skin inflammation.

Application Notes

Optimal dilution of the CCR10 antibody should be determined by the researcher.

Immunogen

Amino acids TEATEQVSWGHYSGDEEDAYSAEPLPELCYKADVQAFSRAFQ from the human protein were used as the immunogen for the CCR10 antibody.

Storage

After reconstitution, the CCR10 antibody can be stored for up to one month at 4oC. For long-term, aliquot and store at -20oC. Avoid repeated freezing and thawing.