


CD10 Antibody / Neprilysin / Cytoplasmic domain [clone RM337] (R20361)

Catalog No.	Formulation	Size
R20361-0.1ML	Antibody in PBS with 50% glycerol, 1% BSA and 0.09% sodium azide	100 ul


Recombinant **RABBIT MONOCLONAL**

Bulk quote request


Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Clonality	Recombinant Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	RM337
Purity	Protein A purified from animal origin-free supernatant
UniProt	P08473
Applications	Immunohistochemistry (FFPE) : 1:200-1:500 Western Blot : 1:1000-1:2000
Limitations	This recombinant CD10 antibody is available for research use only.



Western blot testing of human Raji cell lysate with recombinant CD10 antibody at 1:1:000. Routinely visualized at ~100 kDa.



IHC staining of FFPE human kidney tissue with recombinant CD10 antibody at 1:200.



IHC staining of FFPE human breast cancer tissue with recombinant CD10 antibody at 1:200.

Description

This antibody reacts to the cytoplasmic domain of human CD10.

Application Notes

The stated application concentrations are suggested starting points. Titration of the recombinant CD10 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A peptide corresponding to the N-terminus of human CD10 was used as the immunogen for the recombinant CD10 antibody.

Storage

Store the recombinant CD10 antibody at -20°C.