


## CD123 Antibody / IL3RA (R35737)

Catalog No.	Formulation	Size
R35737-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3563
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This CD123 antibody is available for research use only.



Western blot testing of Jurkat [A] and MOLT4 [B] lysates with CD123 antibody at 2ug/ml. Predicted molecular weight ~43kDa.

### Description

Additional name(s) for this target protein: Interleukin 3 receptor alpha; IL3RA

### Application Notes

Optimal dilution of the CD123 antibody should be determined by the researcher.

### Immunogen

Amino acids RQQYECLHYKTD were used as the immunogen for this CD123 antibody.

## Storage

Aliquot and store the CD123 antibody at -20oC.