


CD19 Antibody / Cytoplasmic domain [clone RM332] (R20356)

Catalog No.	Formulation	Size
R20356-0.1ML	Antibody in PBS with 50% glycerol, 1% BSA and 0.09% sodium azide	100 ul


Recombinant **RABBIT MONOCLONAL**

[Bulk quote request](#)


Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Clonality	Recombinant Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	RM332
Purity	Protein A purified from animal origin-free supernatant
UniProt	P15391
Localization	Cell surface, cytoplasmic
Applications	Immunohistochemistry (FFPE) : 1:500-1:1000 Western Blot : 1:1000-1:2000
Limitations	This recombinant CD19 antibody is available for research use only.




Western blot testing of human Raji cell lysate with CD19 antibody at 1:1000. Expected molecular weight: 60~100 kDa depending on glycosylation level.



IHC staining of FFPE human tonsil tissue with recombinant CD19 antibody at 1:1000.



IHC staining of FFPE human tonsil tissue with recombinant CD19 antibody at 1:1000.



IHC staining of FFPE human colon tissue with recombinant CD19 antibody at 1:1000.

Description

This antibody reacts to the cytoplasmic domain of human CD19.

Application Notes

The stated application concentrations are suggested starting points. Titration of the recombinant CD19 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A peptide corresponding to the C-terminus of human CD19 was used as the immunogen for the recombinant CD19 antibody.

Storage

Store the recombinant CD19 antibody at -20°C.