


CD1a Antibody [clone RM393] (R20409)

Catalog No.	Formulation	Size
R20409-0.1ML	Antibody in PBS with 50% glycerol, 1% BSA and 0.09% sodium azide	100 ul


Recombinant **RABBIT MONOCLONAL**

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Clonality	Recombinant Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	RM393
Purity	Protein A purified from animal origin-free supernatant
UniProt	P06126
Localization	Cell surface, cytoplasmic
Applications	Immunohistochemistry (FFPE) : 1:100-1:200
Limitations	This recombinant CD1a antibody is available for research use only.



IHC staining of FFPE human thymus tissue with recombinant CD1a antibody at 1:100.



IHC staining of FFPE human skin tissue with recombinant CD1a antibody at 1:100.

Description

This antibody reacts to human CD1a.

Application Notes

The stated application concentrations are suggested starting points. Titration of the recombinant CD1a antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A peptide corresponding to Human CD1a was used as the immunogen for the recombinant CD1a antibody.

Storage

Store the recombinant CD1a antibody at -20oC.