


CD25 Antibody / IL2RA [clone AFGC-9] (RQ4974)

Catalog No.	Formulation	Size
RQ4974	Antibody in PBS with 0.02% sodium azide, 50% glycerol and 0.4-0.5mg/ml BSA	100 ul

[Bulk quote request](#)

Availability	1-2 weeks
Species Reactivity	Human
Format	Purified
Clonality	Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	AFGC-9
Purity	Affinity purified
UniProt	P01589
Applications	Immunohistochemistry (FFPE) : 1:50-1:200
Limitations	This CD25 antibody is available for research use only.



IHC staining of FFPE human tonsil tissue with CD25 antibody. HIER: boil tissue sections in pH6, 10mM citrate buffer, for 10-20 min and allow to cool before testing.

Description

Receptor for interleukin-2. The receptor is involved in the regulation of immune tolerance by controlling regulatory T cells (TREGs) activity. [UniProt]

Application Notes

Optimal dilution of the CD25 antibody should be determined by the researcher.

Immunogen

A synthetic peptide specific to human CD25 / IL2RA was used as the immunogen for the CD25 antibody.

Storage

Store the CD25 antibody at -20oC.