


CD27 Antibody (R35484)

Catalog No.	Formulation	Size
R35484-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	939
Localization	Cytoplasmic, membranous
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This CD27 antibody is available for research use only.



Western blot testing of human spleen [A], tonsil [B] and thymus [C] lysates with CD27 antibody at 0.3ug/ml. Predicted molecular weight ~29/45kDa (unmodified/glycosylated)

Description

Additional name(s) for this target protein: T-cell activation antigen CD27, T cell activation antigen S152

Application Notes

Optimal dilution of the CD27 antibody should be determined by the researcher.

Immunogen

Amino acids QTLADFRQLPARTLS were used as the immunogen for this CD27 antibody.

Storage

Aliquot and store the CD27 antibody at -20oC.