


## CD29 Antibody / Integrin beta 1 / ITGB1 (R34833)

| Catalog No.  | Formulation                                                                   | Size   |
|--------------|-------------------------------------------------------------------------------|--------|
| R34833-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

|                    |                                                             |
|--------------------|-------------------------------------------------------------|
| Availability       | 1-3 business days                                           |
| Species Reactivity | Human                                                       |
| Format             | Antigen affinity purified                                   |
| Clonality          | Polyclonal (goat origin)                                    |
| Isotype            | Goat Ig                                                     |
| Purity             | Antigen affinity                                            |
| Gene ID            | 3688                                                        |
| Applications       | Western Blot : 0.3-1ug/ml<br>ELISA (peptide) LOD : 1:128000 |
| Limitations        | This CD29 antibody is available for research use only.      |



Western blot testing of gastrointestinal cancer [A] and gall bladder cancer [B] lysates with CD29 antibody at 0.5ug/ml. Predicted molecular weight: 88/130~150 kDa (unmodified/glycosylated).

### Description

Additional name(s) for this target protein: Integrin beta 1; ITGB1

### Application Notes

Optimal dilution of the CD29 antibody should be determined by the researcher.

### Immunogen

Amino acids NKCPKKDSDFKIRP were used as the immunogen for this CD29 antibody.

## Storage

Aliquot and store the CD29 antibody at -20oC.