


CD30 Antibody (F52777)

Catalog No.	Formulation	Size
F52777-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F52777-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	P28908
Applications	Western Blot : 1:1000
Limitations	This CD30 antibody is available for research use only.



Western blot analysis of lysate from HeLa cell line using CD30 antibody diluted at 1:1000. Predicted molecular weight: 53-120 kDa depending on glycosylation level.

Description

Receptor for TNFSF8/CD30L. May play a role in the regulation of cellular growth and transformation of activated lymphoblasts. Regulates gene expression through activation of NF- kappa-B.

Application Notes

Titration of the CD30 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

This CD30 antibody was produced from a rabbit immunized with a KLH conjugated synthetic peptide between 30-63 amino acids from the N-terminal region of human TNFRSF8/CD30.

Storage

Aliquot the CD30 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.