


Cd4 Antibody (RQ6834)

Catalog No.	Formulation	Size
RQ6834	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Rat
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit IgG
Purity	Antigen affinity purified
Buffer	Lyophilized from 1X PBS with 2% Trehalose
UniProt	P05540
Applications	Immunohistochemistry (FFPE) : 2-5ug/ml Direct ELISA : 0.1-0.5ug/ml
Limitations	This Cd4 antibody is available for research use only.



IHC staining of FFPE rat spleen tissue with Cd4 antibody. HIER: boil tissue sections in pH8 EDTA for 20 min and allow to cool before testing.

Description

The CD4 antigen acts as a coreceptor with the T-cell receptor on the T lymphocyte to recognize antigens displayed by an antigen presenting cell in the context of class II MHC molecules.

Application Notes

Optimal dilution of the Cd4 antibody should be determined by the researcher.

Immunogen

Recombinant rat protein (amino acids K28-A246) was used as the immunogen for the Cd4 antibody.

Storage

After reconstitution, the Cd4 antibody can be stored for up to one month at 4oC. For long-term, aliquot and store at -20oC. Avoid repeated freezing and thawing.