


CD57 Antibody (R34315)

Catalog No.	Formulation	Size
R34315-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Mouse, Pig, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	27087
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This CD57 antibody is available for research use only.



Western blot testing of human frontal cortex lysate with CD57 antibody at 0.1ug/ml with [B] and without [A] immunizing peptide. Predicted/observed molecular weight ~38kDa.

Description

Additional name(s) for this target protein: HNK1; B3GAT1

Application Notes

Optimal dilution of the CD57 antibody should be determined by the researcher.

Immunogen

Amino acids HVETPRNYKLRGD were used as the immunogen for this CD57 antibody.

Storage

Aliquot and store the CD57 antibody at -20oC.